



## Mining Proteomics Data

### Protein expression in *P. falciparum* apicoplast vs. ER

#### Learning objectives

- Explore proteomics data on VEuPathDB
- Perform a quantitative mass spec. evidence search


#### Introduction

The malaria parasite *Plasmodium falciparum* and other *Plasmodium* species contain a non-photosynthetic plastid organelle called the apicoplast that is crucial to the malaria parasite's survival. Due to the algal origin of the apicoplast (which contains its own DNA), many proteins and pathways are not shared by the human host, making it an attractive target for antimalarial drugs.

To identify genes that are unique to the apicoplast and not present in human hosts, we could ask the question: **What genes in *Plasmodium falciparum* have protein products that are present at a higher concentration in the apicoplast compared to the endoplasmic reticulum (ER)?** Furthermore, we can leverage protein subcellular localization data to limit the results to the apicoplast.

Overview of the search strategy: Use **Quantitative Mass Spec. Evidence** search in two steps-

- (1) identify genes that are upregulated in apicoplast vs. ER
- (2) intersect that search with an apicoplast localization search




1. Navigate to [PlasmoDB.org](https://PlasmoDB.org)




2. From the sidebar or header, search or scroll for “quantitative” and click on **Quantitative Mass Spec. Evidence**

3. Run a search to find genes that are upregulated in the apicoplast sample compared to the ER sample.

- Select the experiment called **Apicoplast and ER Proteomes (Quantitative)(Dd2) (Boucher et al)**- it compares the samples we are interested in.
- Choose the fold change (FC) search.
- Configure this search to return all genes that are upregulated by 1.5 fold in the apicoplast sample compared to the ER sample.



How many genes did you get that have upregulated protein expression in apicoplast compared to ER?



4. Limit the results by leveraging available subcellular localization data. PlasmoDB has a data set that returns genes with the transit peptides that mediate protein targeting to the apicoplast.
  - a. Click on the **add step** button and find the **subcellular localization** search
  - b. Make sure Apicoplast localization is selected and click on the Run Step button

The image shows three screenshots of the PlasmoDB search strategy interface. The top-left screenshot shows a search strategy titled "Upregulated in apicoplast vs. ER \*" with Step 1: "PfDd2 Apico ER Prot (fc) 272 Genes". A red circle highlights the "Add a step" button. The top-right screenshot shows the "Add a step to your search strategy" dialog with the "Combine with other Genes" section selected. A red arrow points to the "1 INTERSECT 2" option. The bottom screenshot shows the "Search for Genes by Pfal 3D7 Subcellular Localization" dialog. A red arrow points to the "Apicoplast" option in the "Localization" dropdown, and another red circle highlights the "Run Step" button.

5. How many genes did you identify? Explore your results.

Are you more confident that these genes are apicoplast genes? How would you use the PlasmoDB tools to boost your confidence in these so-called apicoplast genes?

